

Chemically Intelligent Metabolite ID Workflows

Dave Heywood
MS Field Marketing

- Metabolite Profiling or Metabolite ID?
- Mass Spectrometry and Accurate Mass Measurements
- Metabolite ID Workflow
 - Comprehensive data collection
 - Intelligent data interpretation
 - Tools for interpreting structure and elemental composition
 - Qual/Quan, integrating quantitative measurements
- Gratuitous commercial type slides

What's Metabolic Profiling at Waters?



Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

Metabolic Profiling

**METABOLITE ID
AND DRUG
METABOLISM**

**BIOMARKER ANALYSIS
LIPIDOMICS AND METABOLOMICS**

**TRADITIONAL
MEDICINE**




Analytical Challenges in Metabolic Profiling

- **Extracting the maximum amount of information**
 - Complex mixtures with wide dynamic range
 - Full structural characterization
- **Providing wide range of experimental options**
 - New ways of extracting more information
- **Increase productivity**
 - Provide consistently high performance for users of all experience levels
 - Provide ability to expand analytical possibilities in the future



Acquity UPLC The Gold Standard in Liquid Chromatography


Waters
THE SCIENCE OF WHAT'S POSSIBLE.™





Speed



Resolution



||



Sensitivity

Better MS Data


Tandem Quadrupole MS

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



- Highest Sensitivity for targeted analysis
- Accurate quantitative data
- True class specific screening
- Affordable rugged mass spectrometry technology

TQ Schematic




Xevo G2 QTof



SYNAPT™ G2



- High resolution mass spectra
 - Resolution independent of scan speed
 - MS and MS/MS
- High mass measurement accuracy (Accurate Mass)
 - MS and MS/MS
- High sensitivity full scan data
- MS^E
 - Structural analysis for metabolite localization
- Ideal for qualitative and more recently combined qualitative and quantitative workflows



What is Accurate Mass?

- carbon has a mass of 12
- hydrogen has a mass of 1
- oxygen has a mass of 16
- nitrogen has a mass of 14

But this is not strictly “Accurate”

- carbon has a mass of 12.0000
 - hydrogen has a mass of 1.0078
 - oxygen has a mass of 15.9949
 - nitrogen has a mass of 14.0031
-
- It is possible to have combinations of atoms which have the same nominal (or integer) mass but different accurate mass
 - If such compounds can be mass measured with sufficient accuracy it is possible to determine elemental composition

- CO = 27.9949
 - N₂ = 28.0061
 - C₂H₄ = 28.0313
-
- These elemental combinations have the same nominal mass but different accurate masses
 - A nominal mass measurement cannot distinguish these
 - If any compounds differ in their elemental compositions by substitution of any of these elements, then the exact mass measurement will show this

- The accuracy of the measurement is quoted as the difference (error) between the measured mass and the calculated mass
- The accuracy is measured in
 - milliDaltons (1mDa = 0.001 mass units)
 - ppm = parts per million = $\Delta m/m \times 10^6$

Example:


True' mass = 400.0000

Measured mass = 400.0020

Difference = 0.0020 (2 mDa)

$$\text{ppm error} = \frac{0.002}{400} \times 10^6 = 5 \text{ ppm}$$


When Exact Mass Makes a Difference (assuming elements $C_{50}H_{100}N_5O_5S_2$)




Nominal mass measured spectrum (ie quadrupole or ion trap data)

~50ppm tolerance @ m/z
311.0814


37 possible results



When Exact Mass Makes a Difference (assuming elements $C_{50}H_{100}N_5O_5S_2$)




Exact Mass Measured Spectrum
3ppm tolerance @ m/z 311.0814
3 Possible Results




Exact mass for Selectivity Isobaric metabolites of Rabeprazole

Loss of O =
-16 mass shift




Sulphide

[M+H]⁺ 344.1433



Difference = 36 mDa

Loss of CH₄ =
-16 mass shift



Aldehyde

[M+H]⁺ 344.1069

What are the Benefits of Exact Mass?


- Measurement of mass to 4 decimal places
- High confidence in confirming expected compounds
 - Distinguishes them from compounds of similar mass
- Compound identification
 - Prediction of elemental composition
- Patent submission and publication
 - ACS require better than 5ppm mass accuracy for publication


What instruments are used for exact mass?


- Magnetic sector mass spectrometers
 - these have traditionally been used for exact mass measurement
 - required skilled operator to get good results
- Orthogonal time-of-flight (oa-TOF) mass spectrometers
 - routine operation with good mass accuracy (<1ppm)
 - Accurate isotope ratio measurements
- Ion Cyclotron Resonance mass spectrometers (FTICR, Orbitrap)
 - generally expensive and more difficult to operate
 - good mass accuracy (<1ppm)

Advanced Elemental Composition Calculations

Elemental Composition Analysis







Exact Mass Measurement Elemental Composition – Isotopic Fit

molecular mass [Da]	without isotope abundance information					2% isotopic abundance accuracy	5% isotopic abundance accuracy
	10 ppm	5 ppm	3 ppm	1 ppm	0.1 ppm	3 ppm	5 ppm
150	2	1	1	1	1	1	1
200	3	2	2	1	1	1	1
300	24	11	7	2	1	1	6
400	78	37	23	7	1	2	13
500	266	115	64	21	2	3	33
600	505	257	155	50	5	4	36
700	1046	538	321	108	10	10	97
800	1964	973	599	200	20	13	111
900	3447	1712	1045	345	32	18	196

BMC Bioinformatics 2006, 7:234 doi:10.1186/1471-2105-7-234

Metabolite ID and Drug Metabolism



- Pharmaceutical R&D is a long, costly and risky activity. On average it takes 12 years to develop and market a NME

The price of innovation: new estimates of drug development costs

Joseph A. DiMasi^{a,*}, Ronald W. Hansen^b, Henry G. Grabowski^c

^a *Tufts Center for the Study of Drug Development, Tufts University, 192 South Street, Suite 550, Boston, MA 02111, USA*

^b *William E. Simon Graduate School of Business Administration, University of Rochester, Rochester, NY, USA*

^c *Department of Economics, Duke University, Durham, NC, USA*

Received 17 January 2002; received in revised form 24 May 2002; accepted 28 October 2002

Abstract

The research and development costs of 68 randomly selected new drugs were obtained from a survey of 10 pharmaceutical firms. These data were used to estimate the average pre-tax cost of new drug development. The costs of compounds abandoned during testing were linked to the costs of compounds that obtained marketing approval. The estimated average out-of-pocket cost per new drug is US\$ 403 million (2000 dollars). Capitalizing out-of-pocket costs to the point of marketing approval at a real discount rate of 11% yields a total pre-approval cost estimate of US\$ 802 million (2000 dollars). When compared to the results of an earlier study with a similar methodology, total capitalized costs were shown to have increased at an annual rate of 7.4% above general price inflation.

Working with the Pharmaceutical Industry

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

For a decade we have partnered with the pharmaceutical industry to bring meaningful impact

“In terms of efficiency for metabolite identification studies, the accurate mass LC/MS^E approach has provided significant gains...
...typical savings using the approach outlined in this paper have been in the range of 13 hours per molecule. As a result the capacity for conducting preliminary metabolite identification experiments has increased by almost an order of magnitude.”

P. R. TILLER, *et al.*
Rapid Commun. Mass Spectrom.
2008; 22: 1053–1061


“UPLC with QToF and MetaboLynx XS provides an empowering platform for our metabolism scientists...With this complete workflow, we can routinely see more metabolites in a single run, and present a more definitive metabolic pathway picture for our clients in less time.”

DR. DAVID JOHNSON,
DIRECTOR OF DMPK,
MICROCONSTANTS


The 9 biggest pharmaceutical companies worldwide* use Waters Metabolite Identification System Solution

*ranked by global prescription drug sales

Drug Metabolism Workflow to Maximize Productivity



A Typical Metabolite ID Experiment




Full Scan MS
Targeted MS/MS or DDA
Precursor Ion Scan – QQQ
Neutral Loss Scan - QQQ

Analyte minus Control (Exact Mass?)
MDF if Exact Mass data

Combine data from multiple data files (instruments?)


Our Metabolite Identification Workflow

For *in-vitro/in-vivo* samples



Metabolite ID Workflow to Maximize Productivity

Analyze



ACQUITY UPLC

Maximum chromatographic resolution, sensitivity, and speed for MS-based studies



QTOF with MS^E

Exact mass analysis with data-rich information



Chemically Intelligent data processing


Identification of expected & unexpected metabolites
Dealkylation / Mass Defect Filter tools



Tools for structural elucidation

Elecomp^E and MassFragment for comprehensive structural elucidation

Interpret





GOAL = Rapid Confident Metabolite Identifications

Collect MS/MS of Parent Compound

Nefazodone, HLM GSH T90

2ndGSH_1101_010 605 (6.222) Cm (605-(599:602+608:613))


2: TOF MS ES+
417



Confirmation of structures → Structure elucidation step

- The goal is to provide the most likely structure for the fragment ions generated and reduce the bottleneck in the processing and rationalization of structural fragment assignment
- This software algorithm assigns structures by taking fragment ion spectra of the drug or compound and using it to automatically calculate fragments based on a series of novel chemically intelligent algorithms
- This approach is based on systematic bond disconnection for the precursor structure instead of the usual 'rule based approach'

Submission

<input checked="" type="radio"/> Structure															
Product ion(s) (Da)	<table border="1"> <tr><td>55.9926</td><td>5</td></tr> <tr><td>56.0237</td><td>2</td></tr> <tr><td>56.0511</td><td>109</td></tr> <tr><td>69.0095</td><td>4</td></tr> <tr><td>77.0400</td><td>59</td></tr> <tr><td>83.0605</td><td>417</td></tr> <tr><td>83.7390</td><td>1</td></tr> </table>	55.9926	5	56.0237	2	56.0511	109	69.0095	4	77.0400	59	83.0605	417	83.7390	1
55.9926	5														
56.0237	2														
56.0511	109														
69.0095	4														
77.0400	59														
83.0605	417														
83.7390	1														
DBE	<input type="text" value="0"/> to <input type="text" value="50"/>														
Electron count	odd: <input type="radio"/> even: <input type="radio"/> both: <input checked="" type="radio"/>														
Maximum H deficit	<input type="text" value="6"/>														
Fragment number of bonds	one: <input type="radio"/> (fastest) two: <input type="radio"/> three: <input type="radio"/> four: <input checked="" type="radio"/> (fast)														
Scoring method	use SMARTS : <input type="radio"/> use scoring function: <input checked="" type="radio"/>														
Scoring function parameters	phenyl: <input type="text" value="8"/> aromatic: <input type="text" value="6"/> multiple: <input type="text" value="4"/> ring: <input type="text" value="2"/> single: <input type="text" value="1"/>														
General parameters	hetero modifier: <input type="text" value="0.5"/> H-penalty: <input type="text" value="0"/> max score: <input type="text" value="16"/>														
Output order by	mass: <input checked="" type="radio"/> intensity: <input type="radio"/>														

mode: +ve neutral -ve
 n top ions (raw only)
 % Int. cutoff (raw only)
 +/- 0.1 0.01
 Filter no-structure results:

Submit

Systematic Bond Disconnection and Exact Mass - MassFragment

Report

Input





Chiral	
ID (job)	18
Mass (Da)	720.3128
Formula	C ₃₇ H ₄₈ N ₆ O ₅ S ₂
DBE	17

Experiment

Product ion(s) (Da)	171.0973 197.0782 227.2023 268.1506 296.1432 426.1859 507.2465 533.2224 721.3160 +/- 0.01 in positive mode, structure filter on
DBE	0 to 50
Electron count	both
Maximum H deficit	6
Fragment number of bonds	4
Scoring	aromatic: 6, multiple: 4, ring: 2, phenyl: 8, other: 1 H-deficit: 0, hetero modifier: 0.5, max score: 16
Order:	intensity
Plot:	show <input type="radio"/> hide <input checked="" type="radio"/>

Results

Results:

Daughter (Da)	Mass (Da)	Mass error (mDa)	Formula	DBE	ΔFormula	Structure(s), score & H-deficit
274.1545	274.1556	-1.1	C ₁₅ H ₂₀ N ₃ O ₂	7.5	C ₁₀ H ₁₃ N ₂ Cl	  ■ S:0.5 B:1 H:0 <input type="checkbox"/> ■ S:6 B:4 H:0 <input type="checkbox"/>
246.1247	246.1243	+0.4	C ₁₃ H ₁₆ N ₃ O ₂	7.5	C ₁₂ H ₁₇ N ₂ Cl	  ■ S:1 B:1 H:0 <input type="checkbox"/> ■ S:1.5 B:2 H:+1 <input type="checkbox"/>

Nomenclature explanation

A color that corresponds to the score: green (low) -> red (high)


score H-diff

■ S:9 B:2 H:+1 report checkbox

bonds broken

Report






Input:

	ID (job)	38
	Mass (Da)	469.2245
	Formula	C ₂₉ H ₃₂ N ₃ O ₂ Cl
	DBE	12


Experiment:

Product ion(s) (Da)	140.0819 154.0957 180.1124 246.1197 274.1526 275.1581 56.0511 77.0400 83.0605 97.0404 -/- 0.01 in positive mode, structure filter 1
DBE	0 to 50
Electron count	both
Maximum H deficit	6
Fragment number of bonds	4
Scoring	aromatic: 6, multiple: 4, ring: 2, phenyl: 8, other: 1 H-deficit: 0, hetero modifier: 0.5, max score: 16
Order:	mass
Plot:	show <input checked="" type="radio"/> hide <input type="radio"/>
Files:	DMX CSV






Results:

 274.1526 \rightarrow + (+0H) 274.1556 (-3.0 mDa) (S:0.5, B:1) C ₁₅ H ₂₀ N ₃ O ₂ (-C ₁₀ H ₁₃ N ₂ Cl)	 246.1197 \rightarrow + (+0H) 246.1243 (-4.6 mDa) (S:1.0, B:1) C ₁₃ H ₁₆ N ₃ O ₂ (-C ₁₂ H ₁₇ N ₂ Cl)	 180.1124 \rightarrow + (-1H) 180.1137 (-1.3 mDa) (S:1.0, B:2) C ₉ H ₁₄ N ₃ O (-C ₁₆ H ₁₉ N ₂ OCl)
 154.0957 \rightarrow + (+1H) 154.0980 (-2.3 mDa) (S:1.0, B:2) C ₇ H ₁₂ N ₃ O (-C ₁₈ H ₂₁ N ₂ OCl)	 140.0819 \rightarrow + (+1H) 140.0824 (-0.5 mDa) (S:1.0, B:2) C ₆ H ₁₀ N ₃ O (-C ₁₉ H ₂₃ N ₂ OCl)	

Nelazodipe, HLM GSH T90 - 2 n(GSH_1101_010 605 6 222) Cm 605-659602-608 613) - 2. TOF MS ES+




Results:

 274.1526 \rightarrow + (+0H) 274.1556 (-3.0 mDa) (S:0.5, B:1) C ₁₅ H ₂₀ N ₃ O ₂ (-C ₁₀ H ₁₃ N ₂ Cl)	 246.1197 \rightarrow + (+0H) 246.1243 (-4.6 mDa) (S:1.0, B:1) C ₁₃ H ₁₆ N ₃ O ₂ (-C ₁₂ H ₁₇ N ₂ Cl)	 180.1124 \rightarrow + (-1H) 180.1137 (-1.3 mDa) (S:1.0, B:2) C ₉ H ₁₄ N ₃ O (-C ₁₆ H ₁₉ N ₂ OCl)
 154.0957 \rightarrow + (+1H) 154.0980 (-2.3 mDa) (S:1.0, B:2) C ₇ H ₁₂ N ₃ O (-C ₁₈ H ₂₁ N ₂ OCl)	 140.0819 \rightarrow + (+1H) 140.0824 (-0.5 mDa) (S:1.0, B:2) C ₆ H ₁₀ N ₃ O (-C ₁₉ H ₂₃ N ₂ OCl)	


Nefazodone, HLM GSH T90
2ndGSH_1101_010 605 (6.222) Cr

2: TOF MS ES+
417




Acquiring a Comprehensive MS and MS/MS Data Set

Conventional Data Directed LC-MS/MS




Precursor Survey Scan

Conventional Data Directed LC-MS/MS




Product Ion Scan
Precursor Ion Selection
Which one(s)?
Ion Transmission Window +/- 2 Da

Product Ion Spectrum Typically very fast (.1-1 second)



Data Directed Analysis LC-MS/MS




MS² is a SERIAL process


MS² is a BIASED process

MS² is a DISCONTINUOUS process

Alternate Scanning LC-MS (LC-MSE) ...monitor precursors




Alternate Scanning LC-MS (LC-MSE) ...monitor fragments (no pre-selection)




Alternate Scanning LC-MS (LC-MS^E)

...precursor and fragment profiles



Alternate Scanning LC-MS (LC-MS^E)

...time resolved mass measurements




All ions all the time

Bateman *et al*, Anal Chem (2002)
Silva *et al*, MCP (2005)


*Patented technology

Metabolite ID Workflow to Maximize Productivity

Analyze



Interpret



ACQUITY UPLC

Maximum chromatographic resolution, sensitivity, and speed for MS-based studies



QTOF with MS^E

Exact mass analysis with data-rich information



Chemically Intelligent data processing

Identification of expected & unexpected metabolites
Dealkylation / Mass Defect Filter tools




Tools for structural elucidation

Elecomp^E and MassFragment for comprehensive structural elucidation

GOAL = Rapid Confident Metabolite Identifications

- Widely applied to chromatographic data sets by Zhang et al, Bateman et al.
 - Integral part of many LC-based metabolite ID workflows
- MetaboLynx XS chemically intelligent MDF that incorporates novel, structure-based dealkylations

Indinavir



$C_{36}H_{48}N_5O_4$

$MH^+ = 614.$ **3706**

Significance of exact mass filtering

Biotransformation	Nominal mass	Accurate mass	Decimal Place shift
+O	+16	+ 15.9949	- 0.0051
+O ₂	+32	+ 31.9898	- 0.0102
-H ₂	-2	- 2.0157	- 0.0157
-CH ₂	-14	- 14.0157	- 0.0157
-Cl+O	-18	- 17.9662	+0.0338
+C ₂ H ₂ O	+42	+ 42.0106	+0.0106
+SO ₃	+80	+ 79.9568	- 0.0432
+C ₆ H ₈ O ₆	+176	+176.0321	+0.0321
+C ₆ H ₈ O ₇	+192	+192.0270	+0.0270
+C ₂ H ₅ NO ₂ S	+107	+107.0042	+0.0042
+C ₁₀ H ₁₅ N ₃ O ₆ S	+305	+305.0682	+0.0682

Phase 1
Metabolism
< 0.04


Phase 2
Metabolism
< 0.07

*The mass shift may be larger if a compound undergoes O-Dealkylation or N-Dealkylation


Mass defect filtering for full scan MS

Compound A : $m/z = 471.2222$




TIC : No Filter



TIC : Filter 70 mDa




Mass Defect Filtering Filtering out false positives




- It is not a linear relationship with mass!
- Fixed linear MDF difficult to automate because risk of metabolic cleavages
- Depends whether S, Cl or Br present
- We can miss important metabolites if filters are not set-up correctly → false negatives
- The C-Heteroatom tool is key to provide the 'correct MDF's

Dealkylated metabolites → Linear Mass Defect for Verapamil

	ID (job)	155
	Mass (Da)	454.2832
	Formula	C ₂₇ H ₃₈ N ₂ O ₄
	DBE	10

-	none
R_12: -C17H24N2O2	167.0994
R_11: -C16H21NO2	196.1260
R_10: -C11H15NO2	262.1729
R_9: -C10H12O2	291.1994
R_8: -CH2O	425.2726
R_7: -CH2O	425.2726
R_6: -CH2O	425.2726
R_5: -CH2O	425.2726
R_4: -CH2	441.2675
R_3: -CH2	441.2675
R_2: -CH2	441.2675
R_1: -CH2	441.2675
R_0: -CH2	441.2675
none	455.2832

Linear Relationship integer mass vs. decimal places




Dealkylated metabolites → Non-Linear Mass Defect for Glyburide

	ID (job)	156
	Mass (Da)	493.1438
	Formula	C ₂₃ H ₂₈ N ₃ O ₅ SCl
	DBE	11

-	none
R_12: -C16H15N2O4SCl	128.0997
R_11: -C15H21N3O3S	171.0134
R_10: -C15H20N2O3S	186.0243
R_9: -C16H14NO2Cl	207.0725
R_8: -C7H12N2O3S	290.0869
R_7: -C8H6NO2Cl	311.1351
R_6: -C8H5O2Cl	326.1460
R_5: -C7H11NO	369.0597
R_4: -C6H11N	397.0546
R_3: -C6H10	412.0655
R_2: -Cl+H	460.1827
R_1: -CH2O	464.1332
R_0: -CH2	480.1281
none	494.1438


Non-Linear Relationship integer mass vs. decimal places



Linear Fixed MDF vs. Intelligent MDF


Combined Metabolite Peaks (All Found and Unexpected Peaks) [Analyte]

1.01e4




Metabolic cleavages detected by Intelligent MDF

Parent compound fragment ion characterization with Metabolynx, MS^E & MassFragment




Dealkylations and other Biotransformations



Metabolite ID Workflow to Maximize Productivity

Analyze



ACQUITY UPLC

Maximum chromatographic resolution, sensitivity, and speed for MS-based studies



Xevo QTOF with MS^E

Exact mass analysis with data-rich information




Chemically Intelligent data processing

Identification of expected & unexpected metabolites
Dealkylation / Mass Defect Filter tools



Tools for structural elucidation

Elecomp^E and MassFragment for comprehensive structural elucidation



Interpret




GOAL = Rapid Confident Metabolite Identifications

IsoCount Metabolite Localization


Verapamil t=30
SMC_010 1115 (6.255) QT (1); Cm (1114:1117-1121:1123)

2: TOF MS ES+
9.39e3

Metabolite Spectrum



Parent Spectrum



IsoCount - Submission


<input checked="" type="radio"/> Structure	
Biotransform	-90.0468 --- "-C7-H6" -78.9105 --- "-Br+H" -67.9984 --- "-C-F3+H" -67.9222 --- "-Cl2+H2" -61.9156 --- "-Br+H+O"
Metabolite ion(s) (Da)	73.0288 1119 118.0868 748 151.0723 3177 188.0902 2285
Parent ion(s) (Da)	118.0862 1776 150.0663 1109 165.0907 21823 167.9989 2328
Mode	<input checked="" type="radio"/> positive <input type="radio"/> negative
Maximum H deficit	4
Mass Tolerance (mDa)	10
Fragment number of bonds	one: <input type="radio"/> (fastest) two: <input type="radio"/>
Scoring method	use SMARTS : <input checked="" type="radio"/>
General parameters	hetero modifier: 0.5 H-penalty: 0
Result parameters	<input type="radio"/> show all <input checked="" type="radio"/> show results for plausible biotransforms Order fragments by <input checked="" type="radio"/> mass <input type="radio"/> intensity

Submit

Spectral Peak Match Mirror Plot

Metabolite Spectrum

Preserved
Shifted
New
Unmatched



Parent Spectrum

Viewing Quantitative Data

MetaboLynx XS Browser - [18Mar10_AFAMM]

File Edit View Tools Window Help

Plate: 2 Vial: 4

Samples, Time Course

Expected Metabolites - 18Mar10_09AFAMM_MDF_25, Verapamil 1 uM T60, MSe, Resolution, pare

Status	m/z Found	Metabolite Name	Formula	mDa	Time
✓	617.3074	Glucuronide conjugation-CH2 (R_0:-CH2)	C32H44N2O10	0.0	6.75
✓	617.3067	Glucuronide conjugation-CH2 (R_0:-CH2)	C32H44N2O10	-0.7	6.37
✓	604.2974	Deethylation + Glucuronide conjugation	C31H43N2O10	-2.2	6.80
✓	489.2976	Alkenes to dihydrodiol	C27H40N2O6	1.2	9.67
✓	471.2850	Hydroxylation	C27H38N2O5	-0.9	9.06
✓	457.2712	Demethylation + hydroxylation	C26H36N2O5	1.0	6.65
✓	457.2681	Demethylation + hydroxylation	C26H36N2O5	-2.1	6.10
✓	455.2901	Parent	C27H38N2O4	-0.9	8.87
✓	441.2740	Demethylation	C26H36N2O4	-1.3	8.70
✓	427.2608	Deethylation	C25H34N2O4	1.1	6.29
✓	427.2600	Deethylation	C25H34N2O4	0.3	8.00
✓	307.2027	Hydroxylation-C10H12O2 (R_5:-C10H12O2)	C17H26N2O3	0.6	4.31
✓	307.2000	Hydroxylation-C10H12O2 (R_5:-C10H12O2)	C17H26N2O3	-2.1	3.29
✓	293.1831	Demethylation + hydroxylation-C10H12O2 (R_5:-C10H12O2)	C16H24N2O3	-3.4	5.27
✓	291.2064	Parent-C10H12O2 (R_5:-C10H12O2)	C17H26N2O2	-0.8	6.53
✓	277.1899	Demethylation-C10H12O2 (R_5:-C10H12O2)	C16H24N2O2	0.0	4.86

250(Time: 6.75) Combine


1:TOF MS ES+ 3.71e4

Combined Metabolite Peaks (All Found Peaks) [Analyte]


Parent-C10H12O2 (R_5:-C10H12O2)

290
6.53

Accessible Quantitative Analysis



Generates snapshot of metabolism for directing future studies




Alprazolam

#	Name	Sample Text	Type	Std. Conc	RT	Area	Response	Primary Flags	ng/mL	%Dev	
1	1	25Mar10_Alp013	0.01ng/mL alprazolam in plasma	Standard	0.010	1.43	8.129	0.011	bb	0.0127	26.5
2	2	25Mar10_Alp014	0.01ng/mL alprazolam in plasma	Standard	0.010	1.42	10.645	0.014	bb	0.0154	53.7
3	3	25Mar10_Alp015	0.01ng/mL alprazolam in plasma	Standard	0.010	1.43	8.788	0.012	bb	0.0134	34.4
4	4	25Mar10_Alp016	QC 0.02ng/mL alprazolam in plasma	Standard	0.020	1.43	16.806	0.022	bb	0.0222	11.1
5	5	25Mar10_Alp017	QC 0.02ng/mL alprazolam in plasma	Standard	0.020	1.43	17.132	0.023	bb	0.0231	15.4
6	6	25Mar10_Alp018	QC 0.02ng/mL alprazolam in plasma	Standard	0.020	1.43	14.285	0.020	bb	0.0203	1.5
7	7	25Mar10_Alp019	0.05ng/mL alprazolam in plasma	QC	0.050	1.43	35.487	0.050	bb	0.0467	-6.6
8	8	25Mar10_Alp020	0.05ng/mL alprazolam in plasma	QC	0.050	1.43	36.876	0.051	bb	0.0473	-5.5
9	9	25Mar10_Alp021	0.05ng/mL alprazolam in plasma	QC	0.050	1.42	35.192	0.050	bb	0.0463	-7.5
10	10	25Mar10_Alp022	0.1ng/mL alprazolam in plasma	Standard	0.100	1.43	74.748	0.106	bb	0.0946	-5.4
11	11	25Mar10_Alp023	0.1ng/mL alprazolam in plasma	Standard	0.100	1.43	72.988	0.104	bb	0.0930	-7.0
12	12	25Mar10_Alp024	0.1ng/mL alprazolam in plasma	Standard	0.100	1.43	66.447	0.099	bb	0.0887	-11.3
13	13	25Mar10_Alp025	QC 0.2ng/mL alprazolam in plasma	Standard	0.200	1.43	131.644	0.195	bb	0.1724	-13.8
14	14	25Mar10_Alp026	QC 0.2ng/mL alprazolam in plasma	Standard	0.200	1.43	126.307	0.190	bb	0.1674	-16.3
15	15	25Mar10_Alp027	QC 0.2ng/mL alprazolam in plasma	Standard	0.200	1.43	132.971	0.198	bb	0.1750	-12.5
16	16	25Mar10_Alp028	0.5ng/mL alprazolam in plasma	QC	0.500	1.43	322.980	0.500	bb	0.4362	-12.8
17	17	25Mar10_Alp029	0.5ng/mL alprazolam in plasma	QC	0.500	1.43	323.932	0.497	bb	0.4340	-13.2
18	18	25Mar10_Alp030	0.5ng/mL alprazolam in plasma	QC	0.500	1.43	317.463	0.511	bb	0.4460	-10.8
19	19	25Mar10_Alp031	1.0ng/mL alprazolam in plasma	Standard	1.000	1.43	692.117	1.096	bb	0.9522	-4.8
20	20	25Mar10_Alp032	1.0ng/mL alprazolam in plasma	Standard	1.000	1.43	705.248	1.096	bb	0.9521	-4.8
21	21	25Mar10_Alp033	1.0ng/mL alprazolam in plasma	Standard	1.000	1.43	675.462	1.082	bb	0.9401	-6.0
22	22	25Mar10_Alp034	QC 2.0ng/mL alprazolam in plasma	Standard	2.000	1.43	1262.481	2.134	bb	1.8510	-7.5
23	23	25Mar10_Alp035	QC 2.0ng/mL alprazolam in plasma	Standard	2.000	1.43	1252.423	2.064	bb	1.7907	-10.5
24	24	25Mar10_Alp036	QC 2.0ng/mL alprazolam in plasma	Standard	2.000	1.43	1209.636	2.112	bb	1.8320	-8.4
25	25	25Mar10_Alp037	5ng/mL alprazolam in plasma	QC	5.000	1.43	3194.469	5.423	bb	4.6994	-6.0
26	26	25Mar10_Alp038	5ng/mL alprazolam in plasma	QC	5.000	1.43	3190.724	5.413	bb	4.6907	-6.2
27	27	25Mar10_Alp039	5ng/mL alprazolam in plasma	QC	5.000	1.43	3150.381	5.421	bb	4.6980	-6.0
28	28	25Mar10_Alp040	10ng/mL alprazolam in plasma	Standard	10.000	1.43	6062.970	10.839	bb	9.3891	-6.1
29	29	25Mar10_Alp041	10ng/mL alprazolam in plasma	Standard	10.000	1.43	6020.000	10.916	bb	9.4557	-5.4
30	30	25Mar10_Alp042	10ng/mL alprazolam in plasma	Standard	10.000	1.43	5743.824	10.956	bb	9.4905	-5.1
31	31	25Mar10_Alp043	QC 20ng/mL alprazolam in plasma	Standard	20.000	1.43	11536.141	21.841	bb	18.9165	-5.4
32	32	25Mar10_Alp044	QC 20ng/mL alprazolam in plasma	Standard	20.000	1.43	11535.521	22.309	bb	19.3223	-3.4
33	33	25Mar10_Alp045	QC 20ng/mL alprazolam in plasma	Standard	20.000	1.43	11409.267	21.878	bb	18.9486	-5.3
34	34	25Mar10_Alp046	50ng/mL alprazolam in plasma	Standard	50.000	1.43	25764.422	55.504	bb	48.0681	-3.9
35	35	25Mar10_Alp047	50ng/mL alprazolam in plasma	Standard	50.000	1.43	25481.125	55.548	bb	48.1059	-3.8
36	36	25Mar10_Alp048	50ng/mL alprazolam in plasma	Standard	50.000	1.43	25182.746	57.220	bb	49.5536	-0.9
37	37	25Mar10_Alp049	100ng/mL alprazolam in plasma	QC	100.000	1.43	44893.945	121.597	bb	105.3020	5.3
38	38	25Mar10_Alp050	100ng/mL alprazolam in plasma	QC	100.000	1.43	44471.485	119.928	bb	103.8566	3.9
39	39	25Mar10_Alp051	100ng/mL alprazolam in plasma	QC	100.000	1.43	44099.945	122.260	bb	105.8763	5.9
40	40	25Mar10_Alp052	QC 200ng/mL alprazolam in plasma	Standard	200.000	1.43	85813.773	227.101	bb	196.6647	-1.7
41	41	25Mar10_Alp053	QC 200ng/mL alprazolam in plasma	Standard	200.000	1.43	85639.305	236.276	bb	204.6101	2.3
42	42	25Mar10_Alp054	QC 200ng/mL alprazolam in plasma	Standard	200.000	1.43	84212.516	240.477	bb	208.2485	4.1

Alprazolam MS^E Data Quantitation Curve

Compound name: Alprazolam
Correlation coefficient: $r = 0.999436$, $r^2 = 0.998873$
Calibration curve: $1.15478 * x + -0.00374489$
Response type: Internal Std (Ref 1), Area * (IS Conc. / IS Area)
Curve type: Linear, Origin: Exclude, Weighting: $1/x$, Axis trans: None



Additional Capabilities

Application System Solutions

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



MetaboLynx™ XS MarkerLynx™ XS BiopharmaLynx™ i-FIT™ ChromaLynx™
TargetLynx™ MassFragment™ OpenLynx™ ProteinLynx Global SERVER™

'Game-Changing' ...accesses the widest range of compounds & applications

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



Multimode source
ESI – Electrospray Ionization
APGi – Atmospheric Pressure
Chemical Ionization
ESCI® – Dual ESI and APGi



Dual mode source
APPI – Atmospheric
Pressure Photo
Ionization
APCi – Atmospheric
Pressure Chemical
Ionization



*TRIZAIC™ Source
with nanoTile Technology.
Plug & Play nanoFlow*



nanoFlow™ ESI

MALDI nanoFlow ESI APCi ESCi APPI TRIZAIC ASAP APGC



*MALDI – Matrix Assisted
Laser Desorption Ionization*



*ASAP – Atmospheric Pressure
Solids Analysis Probe*




*APGC – Atmospheric Pressure
Gas Chromatography*

- ASAP
 - Atmospheric pressure solids analysis probe
- Direct sample analysis
 - Fast analysis
 - No sample prep
 - No chromatography
 - Solids and liquids



Waters ASAP Metabolites in Urine

Urine sample from patient after dosing with paracetamol (1000mg) and dihydrocodeine (30mg)



Sample capillary loaded with 1µL neat urine

New Universal source options

- APGC
 - Atmospheric pressure GC interface
- Extend compound coverage
 - LC & GC on one instrument
 - Very simple to exchange ion sources
 - Clean APCI type spectra




- Xevo G2 QTOF or Synapt G2 collect sensitive accurate mass data on ALL precursors and products eliminating extra analysis on the same sample
- MetaboLynx XS offers proven, intelligent (structure driven) interpretation of data in an Metabolite ID workflow
- Advanced structural interpretation and advanced elemental composition calculations allow for conclusive metabolite confirmation


**Rapid, Confident, Metabolite
Identifications**

Innovation. Productivity. Effective Decision Making

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



Isolate



Identify

Localize

Report

